

Polyisoprene Under Powder free surgical gloves

Date: 31st March, 2023

Medispo® Polyisoprene Under

Material: Polyisoprene
 Dimensions
 Minimum Values
 Meet EN455:2

Size	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0
Length(mm)	290	290	290	290	290	290	290	290
Width(mm)	72±4	77±5	83±5	89±5	95±5	102±6	108±6	114±6

Thickness(single wall)
 Textures
 Average values

Middle finger(mm)	0.21±0.03
Palm(mm)	0.20±0.03
Cuff(mm)	0.15±0.02

Physical Properties
 The median Values

	Before Ageing	After Ageing
Force at Break(N)	≥ 9	≥ 9

AQL
 Colour
 Shape
 Cuff
 External Surface
 Powder amount
 Packaging
 CE-Classification
 Compliances
 Storage Instructions

Final Release, AQL0.65(Water leak test)
 Blue
 Anatomic Shape (L/R) with curved fingers
 Beaded
 Micro-textured
 ≤2mg/dm²
 1pair/ Film pouch, 50pairs/box, 4boxes/carton
 Class IIa
 CE mark, ISO 13485, EN 455:1:2:3.
 Keep out of direct sunlight, store in a cool and dry place
 Keep away from sources of ozone or ignition

Shelf Life
 Sterile

3 Years
 E-Beam

Zhou dingshan
 QC manager



TO BE THE GLOBAL INTELLIGENT LATEX CENTER



Guilin HBM Health Protections, Inc.
 No.1-2, Shuijing East Road, Economic and
 and Technological Development Area, Guilin, China
 Tel:86-773 255 0119 Fax:86-773 255 8110
 Email: guilinhbm@hbmchina.com
 website: www.hbmchina.com